Saturno 300

every patient is a STAT


CRONY[®]
INSTRUMENTS srl

F2 → OPERATION

F3 → ARCHIVE

F4 → METHODS

F5 → PROFILES

F6 → STANDARDS

F7 → CONTROLS

F8 → MAINTENANCE

F9 → OPTIONS

MAIN WINDOW

- **OPERATION** Will guide the operator step by step through all operations.
- **ARCHIVE** Can store thousands of patient data for future consultation online and to print patient reports.
- **METHODS** To introduce an unlimited number of user programmable methods.
- **PROFILES** To program an unlimited number of profiles or organ panels.
- **STANDARDS** To program an unlimited number of standards, calibrators and multi-standard curves.
- **CONTROLS** To introduce an unlimited number of controls for the Quality Control program.
- **MAINTENANCE** On-Line control of all operating steps and procedures in the system. On-Line trouble shooting.
- **OPTIONS** Special personalization of ones analyzer.

METHODOLOGIES



- Open system, unlimited number of methods are all user programmable.
- Only 300 µl of reagent per test.
- Automatic sample pre-dilution if necessary.
- Two reagent chemistries and immunoassays.

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

REAGENT CONTROL

- 30 On-Line reagents individually cooled or left at room temperature.
- On-Line reagent inventory and display of the number of tests available in each of the reagent containers.
- On-Line reagent stability check.
- Bar code reader for reagent identification (optional).

PATIENT ENTRY

- 40 On-Line primary tubes or sample cups with continuous sample loading and programming while analyzer is in operation.
- STATS will be processed within 6 minutes, without disrupting routine testing.
- Bidirectional interface enables to receive and transmit data from and to outside computer.
- Bar code reader for sample identification (optional).

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

QUALITY CONTROL GRAPH

A 60 days ON-LINE quality control with its actual values.

TYPICAL KINETIC MEASUREMENT

Graph that can be displayed for every sample.

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SATURNO 300 is a fast random access analyzer for Clinical Chemistry and Immunoassay, where

EVERY PATIENT IS A STAT.

SATURNO 300 is a 24 hours on duty to run your daily routine, a single test, a profile, an electrolyte, an organ panel or a stat with the same speed and efficiency.

State of the art technology with sophisticated user friendly software offers a new dimension in versatility, speed and reliability. It enables each user to personalize the analyzer to fit his particular needs and become an indispensable collaborator in any laboratory.

300 tests per hour plus 160 ISE electrolytes at the same time.

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
OPERATIONAL CYCLE

- Reagent is transferred into a reaction cuvette and brought to temperature prior to the addition of the sample.
- Cuvette with its reagent is immediately measured in the built-in photometer to obtain a Reagent/Cuvette Blank.
- The sample is transferred and mixed into the previously measured reagent/cuvette to start reaction.
- The reaction cuvettes are immersed in a liquid bath at 37° C and are automatically washed for a non stop operation.
- Sample is transferred to ISE module, if required.
- Automatic print-out of patient report.

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Saturno 300

HIGHLIGHTS

- Continuous patient programming while analyzer is in operation.
- 30 On-Line reagents individually cooled or left at room temperature.
- Cost effective operation - no disposables
- Walk away system - the reaction cuvettes are automatically washed and dried to assure a non stop, continuous operation.
- Liquid sensor probes for sample preparation and reagent dispensing with built-in mixer.
- Check and control of reagent inventory and stability, flagged when reagent is missing or expired.
- Automatic repetition of critical or out of range results.
- I.S.E. module for simultaneous determination of Sodium, Potassium and Chloride (optional).
- Identification of sample and reagent with a built-in bar code reader (optional).
- Help on-line.

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TECHNICAL SPECIFICATION**SPECIFICHE TECNICHE**

OPERATING	True Random Access
ASSAY METHODS	Endo point, Fixed Time, Kinetic, Bichromatic, Differential, Immunoassays using Multistandard, Log/logit, or Cubic calibrations.
30 REAGENTS ON-LINE	Individually cooled or at Room Temp.
CHEMISTRY MENU	Open System with unlimited number of Methods, Profiles, Standards, Calibrators and Controls - user programmable.
TROUGHPUT	300 Chemistry tests plus 160 ISE's/Hour ($\pm 5\%$)
RE-RUN OF TESTS	Automatic or upon Request
SAMPLE TRAY	40 On-Line samples on 4 racks in primary tubes or plastic cups that can be continuously added for up to 240 patients at a time. Each patient sample has a corresponding plastic cup for automatic predilution if necessary. 20 extra places are reserved for standards, calibrators, control and STATS.
SAMPLE VOLUME	3 μ l to 50 μ l (with 1 μ l increments)
REAGENT VOLUME	R ₁ from 3 to 500 μ l (with 1 μ l increments) R ₂ from 3 to 250 μ l (with 1 μ l increments)
REAGENT CONSUMPTION	Only 300 μ l per Test
PHOTOMETER	Analytical Plate with 39 Quartz Cuvettes Immersed in a 37°C liquid bath with a Built-in Multi-filter Photometer; Measures the Samples at 8 Wavelength continuously during their incubation. A built-in automatic cuvette washing system assures a non-stop operation.
LIGHT SOURCE	Double halogen lamp - 20 Watt
DETECTOR	Silicon photodiode
CUVETTE VOLUME	600 μ l
OPTICAL PATH	7 mm
PHOTOMETRIC LINEARITY	From -500 to 2.800 Abs $\pm 1\%$
PHOTOMETRIC SENSIVITY	0.0005 Abs

CHARACTERISTICS:

- Industrial Pentium Computer 233 Mhz, 64 Mb RAM, Hard Disk > 20 Gbyte, Floppy disk 1.44 Mb, double serial port.
- Bidirectional Interface to EDP systems.
- Innovating Graphic Software in Window on a large Colored Screen.
- Liquid sensor Probes for Sample Preparation and Reagent Dispensing with built-in Mixer.
- Check and Control of Reagent Inventory and Stability with Signaling when reagent is missing or expired.
- Power Supply - 220/117 Vac 50/60 Hz - 500 W.
- Dimensions - 103x67x45 cm.
- Weight - 104 Kg.

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CE Technical Specifications are subject to change without prior notice

CRONY®
INSTRUMENTS s.r.l.

Via Tiburtina, 661 - 00159 Roma - Italy
Tel. +39 06 43 96 432 - 06 43 81 163 - Fax +39 06 43 90 487
E-mail: mail@crony.com - Internet: www.crony.com